



Detail: PG 8583 CD Laboratory Glassware Washer
Standard Depth

Miele Professional - 9 Independence Way, Princeton, NJ 08540

Detail #:	
Drawn By: Miele Professional	
Date: 9/8/16	Scale: 1"=1'-0"



- EL: Electrical Connection
- KW: Cold Water Connection
- DK: Cold Water Connection/Steam Condenser
- NW: Network and Printer Connection
- WW: Hot Water Connection
- AW: Drain Water
- VE: Demineralized Water Connection
- PA: Equipotential Bonding



Detail: **PG 8583 CD Laboratory Glassware Washer**
- Connections

Miele Professional - 9 Independence Way, Princeton, NJ 08540

Detail #:	
Drawn By: Miele Professional	
Date: 9/8/16	Scale: NTS