


TRUMPF MEDICAL - A HILL-ROM COMPANY

1046 LEGRAND BLVD  
CHARLESTON, SC 29492  
TOLL FREE: 1-888-474-9359  
PHONE: (843) 534-0606  
FAX: (843) 543-0206  
[www.trumpfmedical.com/us](http://www.trumpfmedical.com/us)

SELECT DESIRED MOUNT:

STANDARD

MOTORIZED BASE AND WIRELESS INTERFACE



**SPECIFICATION**

**MAXIMUM LOAD:** 1,000 LBS (450 KG)

**HEIGHT ADJUSTMENT:** 23 IN – 45 IN ( $\pm 0.4$  IN) (59 CM – 114 CM ( $\pm 1$  CM))

**LIFT:** 22 IN  $\pm 0.4$  IN (55 CM  $\pm 1$  CM)

**TRENDELENBURG / REVERSE TRENDELENBURG:**  $\pm 45^\circ$

**TILT:**  $\pm 30^\circ$

**LONGITUDINAL SLIDE:** 18 IN  $\pm 0.4$  IN (46 CM  $\pm 1$  CM)

**ADJUSTMENT RANGE OF LEG SECTION JOINTS :**  $+ 90^\circ / - 105^\circ$


**ADJUSTMENT OF LOWER BACK SECTION JOINTS:**  $+ 90^\circ / - 45^\circ$

**MANUFACTURER NOTES:**

**COPYRIGHT:** [[2018]] HILL-ROM SERVICES, INC. ALL RIGHTS RESERVED

**NOTES:**

1. INSTALLATION TO BE COMPLETED IN ACCORDANCE WITH MANUFACTURER'S SPECIFICATIONS.
2. DO NOT SCALE DRAWING.
3. THIS DRAWING IS INTENDED FOR USE BY ARCHITECTS, ENGINEERS, CONTRACTORS, CONSULTANTS AND DESIGN PROFESSIONALS FOR PLANNING PURPOSES ONLY. THIS DRAWING MAY NOT BE USED FOR CONSTRUCTION.
4. ALL INFORMATION CONTAINED HEREIN WAS CURRENT AT THE TIME OF DEVELOPMENT BUT MUST BE REVIEWED AND APPROVED BY THE PRODUCT MANUFACTURER TO BE CONSIDERED ACCURATE.
5. CONTRACTOR'S NOTE: FOR PRODUCT AND COMPANY INFORMATION VISIT [www.CADdetails.com/info](http://www.CADdetails.com/info) AND ENTER REFERENCE NUMBER 5249-008



**OR SYSTEM TABLES**

TRUSYSTEM® 7000 OR TABLE